

MALT1 Rabbit mAb

货号: B30819

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	91kDa/92kDa
实测分子量	92kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human MALT1 .
序列	PGIKITVNPESKAVLAGQFVKLCCRATGHPFVQYQWFKMNKEIPNGNTSELIFNAVHVKDAGFYVCRVNNNFTFEFSQWS QLDVCDIPESFQRSVDGVSESKLQICVEPTSQLMPGSTLVLCQVAVGSPIPHYQWFKNELPLTHETKKLYMVPYVDLEHQ GTYWCHVYNDRDSQDSKKVEIIDELENNLGHDPNKEQTTDQPLAKDKVALLIGNMNYREHPKLKAPLVDVYELTNLLRQLD FKVVSLLDLTEYEMRNAVDEFLLLLDKGVYGLLYAGHGYENFGNSFMVPVDAPNPYRSENCCLCVQNILKLMQEKETGLNV FLL

靶点信息

研究背景	This gene has been found to be recurrently rearranged in chromosomal translocation with two other genes - baculoviral IAP repeat-containing protein 3 (also known as apoptosis inhibitor 2) and immunoglobulin heavy chain locus - in mucosa-associated lymphoid tissue lymphomas. The protein encoded by this gene may play a role in NF-kappaB activation. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene.
基因ID	10892
基因名	MALT1
Swiss	Q9UDY8
别名	MALT1;IMD12;MLT;MLT1;PCASP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn