

# Stathmin 1 Rabbit mAb

货号: B30896

产品信息

|       |                                                                                                                                       |
|-------|---------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                                       |
| 宿主    | Rabbit                                                                                                                                |
| 克隆性   | Monoclonal                                                                                                                            |
| 预测反应  |                                                                                                                                       |
| 应用    | WB IHC ICC IP FC                                                                                                                      |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50<br><b>FC:</b> 1:20 - 1:50 |
| 理论分子量 | 17kDa/19kDa                                                                                                                           |
| 实测分子量 | 18kDa                                                                                                                                 |
| 形式    | Liquid                                                                                                                                |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                                         |
| 偶联物   | Unconjugated                                                                                                                          |
| 阳性对照  | Jurkat,C6                                                                                                                             |
| 细胞定位  | Cytoplasm,cytoskeleton                                                                                                                |
| 纯化    | Affinity purification                                                                                                                 |

抗原信息


|      |                                                               |
|------|---------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human Stathmin 1. |
| 序列   | ILSPR                                                         |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene belongs to the stathmin family of genes. It encodes a ubiquitous cytosolic phosphoprotein proposed to function as an intracellular relay integrating regulatory signals of the cellular environment. The encoded protein is involved in the regulation of the microtubule filament system by destabilizing microtubules. It prevents assembly and promotes disassembly of microtubules. Multiple transcript variants encoding different isoforms have been found for this gene. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                                           |
|-------|-----------------------------------------------------------|
| 基因ID  | 3925                                                      |
| 基因名   | STMN1                                                     |
| Swiss | P16949                                                    |
| 别名    | STMN1;C1orf215;LAP18;Lag;OP18;PP17;PP19;PR22;SMN;stathmin |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)