

ROC1 Rabbit mAb

货号: B30923

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	12kDa
实测分子量	12kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ROC1 .
序列	MAAAMDVDTPSGTNSGAGKKRFEVKKWNALWAWDIVVDNCAICRNHIMDLCECQANQASATSEECTVAWGVCNH AFHFCISRWLKTRQVCPLDNREWEFQKYGH

靶点信息

研究背景	This locus encodes a RING finger-like domain-containing protein. The encoded protein interacts with cullin proteins and likely plays a role in ubiquitination processes necessary for cell cycle progression. This protein may also affect protein turnover. Related pseudogenes exist on chromosomes 2 and 5.
------	--

基因ID	9978
基因名	RBX1
Swiss	P62877
别名	RBX1;BA554C12.1;RNF75;ROC1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn