

NAT1 Rabbit mAb

货号: B30946

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	33kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,SH-SY5Y
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NAT1.
序列	MDIEAYLERIGYKKS RNKLDLETLDILQHQIRAVPFENLNIHCGDAMD LGLEAIFDQVVRNRN RGGWCLQVNHLLYWALTTI GFETTMLGGYVYSTPAKK

靶点信息

研究背景	This gene is one of two arylamine N-acetyltransferase (NAT) genes in the human genome, and is orthologous to the mouse and rat Nat2 genes. The enzyme encoded by this gene catalyzes the transfer of an acetyl group from acetyl-CoA to various arylamine and hydrazine substrates. This enzyme helps metabolize drugs and other xenobiotics, and functions in folate catabolism. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	9
基因名	NAT1
Swiss	P18440
别名	NAT1;AAC1;MNAT;NAT-1;NATI

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn