

# FDFT1 Rabbit mAb

货号: B30951

## 产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Monoclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC IF/ICC IP                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50 |
| 理论分子量 | 48kDa                                                                                                          |
| 实测分子量 | 48kDa                                                                                                          |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                  |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | HeLa,A-549,U-251MG,Mouse kidney,Rat liver                                                                      |
| 细胞定位  | Endoplasmic reticulum membrane,Multi-pass membrane protein                                                     |
| 纯化    | Affinity purification                                                                                          |

## 抗原信息


|      |                                                                                                          |
|------|----------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human FDFT1.                                                 |
| 序列   | MGIGMAEFLDKHVTSEQEWDKYCHYVAGLVGIGLSRLFSASEFEDPLVGEDTERANSMGLFLQKTNIIRDYLEDQGGRE<br>FWPQEVWSRYVKKLGDFAKPE |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                  |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a membrane-associated enzyme located at a branch point in the mevalonate pathway . The encoded protein is the first specific enzyme in cholesterol biosynthesis, catalyzing the dimerization of two molecules of farnesyl diphosphate in a two-step reaction to form squalene. [provided by RefSeq, Ju l 2008] |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                     |
|-------|---------------------|
| 基因ID  | 2222                |
| 基因名   | FDFT1               |
| Swiss | P37268              |
| 别名    | DGPT; ERG9; SQS; SS |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)