

MAD3 Rabbit mAb

货号: B31000

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,DU145,Jurkat
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human MAD3.
序列	RRARMHIQKLEDQEQRARQLKERLRSKQQSLQRQLEQLRGLAGAAERERLRADSLDSSGLSSERSDSDQEELEVDSVSLVFGGEAELLRGF

靶点信息

研究背景	This gene encodes a member of the Myc superfamily of basic helix-loop-helix leucine zipper transcriptional regulators. The encoded protein forms a heterodimer with the cofactor MAX which binds specific E-box DNA motifs in the promoters of target genes and regulates their transcription. Disruption of the MAX-MXD3 complex is associated with uncontrolled cell proliferation and tumorigenesis. Transcript variants of this gene encoding different isoforms have been described.[provided by RefSeq, Dec 2008]
------	---

基因ID	83463
基因名	MXD3
Swiss	Q9BW11
别名	BHLHC13; MAD3; MYX

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn