

STXBP1 Rabbit mAb

货号: B31001

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	67kDa/68kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,Mouse lung,Mouse brain
细胞定位	Cytoplasm,Membrane,Peripheral membrane protein,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human STXBP1.
序列	HIAEVSQEVTRSLKDFSSSKRMNTGEKTTMRDLSQMLKKMPQYQKELSKYSTHLHLAEDCMKHYQGTVDKLCRVEQDLA MGTDAEGEKIKDPMRAIVPILLDANVSTYDKIRIILLYIFLKNGITTEENLNKLIQHAQIPPEDSEIITNMAHLGVPIVTDSTLRRRS KPERKERISEQTYQLSRWTPIIKDIMEDTIEDKLDTKHYPYISTRSSASFSTTAVSARYGHWKKNKAPGEYRSGPRLIIFILGG VSLNEMRCAYEVTQANGKWEVLIGSTHILTPQKLLDTLKKLNKTDEEISS

靶点信息

研究背景	This gene encodes a syntaxin-binding protein. The encoded protein appears to play a role in release of ne urotransmitters via regulation of syntaxin, a transmembrane attachment protein receptor. Mutations in th is gene have been associated with infantile epileptic encephalopathy-4. Alternatively spliced transcript va riants have been described.
------	--

基因ID	6812
基因名	STXBP1
Swiss	P61764
别名	STXBP1;MUNC18-1;NSEC1;P67;RBSEC1;UNC18

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn