

Glutaminase Rabbit mAb

货号: B31013

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	17kDa/65kDa/73kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,NIH/3T3,Mouse brain,Rat kidney
细胞定位	Cytoplasm,Mitochondrion,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Glutaminase.
序列	KVNPFPKDRWNNTPMDEALHFGHHDVFKILQEYQVQYTPQGSDNGKENQTVHKNLDGLL

靶点信息

研究背景	This gene encodes the K-type mitochondrial glutaminase. The encoded protein is an phosphate-activated amidohydrolase that catalyzes the hydrolysis of glutamine to glutamate and ammonia. This protein is primarily expressed in the brain and kidney plays an essential role in generating energy for metabolism, synthesizing the brain neurotransmitter glutamate and maintaining acid-base balance in the kidney. Alternative splicing results in multiple transcript variants.
基因ID	2744

基因名	GLS
Swiss	O94925
别名	GLS;AAD20;GAC;GAM;GLS1;KGA

产品验证



Western blot analysis of Glutaminase expressed in Mouse brain,Rat kidney,U-251MG using Glutaminase Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn