

PRDX1 Rabbit mAb

货号: B31021

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	22kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PRDX1 .
序列	MSSGNAKIGH PAPNFKATAVMPDGQFKDISLSDYKGKYVVFFFYPLDFTFVCPT EIIAFSDRAEEFKKLNCQVIGASVDSHF CHLAWVNTPKKQGGLGPMNIPLVSDPKRTIAQDYGV LKADEGISFRGLFIIDDKGILRQITVNDLPVGRSVDETLRLVQAFQ FTDKHGEVCPAGWKPGSDTIKPDVQKSKEYFSKQK

靶点信息

研究背景	This gene encodes a member of the peroxiredoxin family of antioxidant enzymes, which reduce hydrogen peroxide and alkyl hydroperoxides. The encoded protein may play an antioxidant protective role in cells, and may contribute to the antiviral activity of CD8(+) T-cells. This protein may have a proliferative effect and play a role in cancer development or progression. Four transcript variants encoding the same protein have been identified for this gene.
------	---

基因ID	5052
基因名	PRDX1
Swiss	Q06830
别名	PRDX1;MSP23;NKEF-A;NKEFA;PAG;PAGA;PAGB;PRX1;PRXI;TDPX2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn