

RANBP3 Rabbit mAb

货号: B31024

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	53kDa/59kDa/60kDa
实测分子量	70,90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,LO2,HT-29,SW620,U-251MG
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human RANBP3.
序列	MADLANEEKPAIAPPVFVFQKDKGQKRSAGGSSPEGGEDSDREDGNYCPPVKRERTSSLTQFPPSQSEERSSGFRCLKPPT LIHGQAPSAGLPSQKPKEQQRSVLRPAVLQAPQPKALSQTVPSSGTNGVSLPADCTGAVPAASPDTAAWRSPSEAADEV CALEEKEPQKNESSNASEEEEACEKKDPATQQAFVFGQNLDRVKLINESVDEADMENAGHP SADTPTATNYFLQYISSSLE NSTNSADASS

靶点信息

研究背景	This gene encodes a protein with a RanBD1 domain that is found in both the nucleus and cytoplasm. This protein plays a role in nuclear export as part of a heteromeric complex. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
------	--

基因ID	8498
基因名	RANBP3
Swiss	Q9H6Z4
别名	RANBP3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn