

UBE4A Rabbit mAb

货号: B31026

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	123kDa
实测分子量	123kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HUVEC,Mouse brain,Mouse kidney,Rat liver,Rat brain
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human UBE4A.
序列	MTDQENNNNNISSNPFAALFGSLADAKQFAAIQKEQLKQSQDELPA SPDDSDNSVSESLDEFDYSVAEISRSFRSQQEICEQLNINHMIQRIFLITLDNSD

靶点信息

研究背景	This gene encodes a member of the U-box ubiquitin ligase family. The encoded protein is involved in mult iubiquitin chain assembly and plays a critical role in chromosome condensation and separation through th e polyubiquitination of securin. Autoantibodies against the encoded protein may be markers for scleroder ma and Crohn's disease. A pseudogene of this gene is located on the long arm of chromosome 3. Alternat ively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	--

基因ID	9354
基因名	UBE4A
Swiss	Q14139
别名	UBE4A; E4; UBOX2; UFD2; ubiquitination factor E4A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn