

# MERTK Rabbit mAb

货号: B31044

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC IP                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50               |
| 理论分子量 | 188kDa                                                                                        |
| 实测分子量 | 188kDa                                                                                        |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | HepG2                                                                                         |
| 细胞定位  | Membrane,Single-pass type I membrane protein                                                  |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human MERTK.                                               |
| 序列   | MGPAPLPLLLGLFLPALWRRRAITEAREEAKPYPLFGPFPGLQTDHTPLLSLPHASGYQPALMFSPQTQGRPHTGNVAIPQ<br>VTSVESKLPPLAFKHTV |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP). [provided by RefSeq, Jul 2008] |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                 |
|-------|---------------------------------|
| 基因ID  | 10461                           |
| 基因名   | MERTK                           |
| Swiss | Q12866                          |
| 别名    | MER; RP38; Tyro12; c-Eyk; c-mer |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)