

hnRNP U Rabbit mAb

货号: B31068

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	88kDa/90kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A431,Jurkat,OVCAR3
细胞定位	Cell surface,Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human hnRNP U.
序列	PKDEDYKQRTQKKAIEVEGKDLPEHAVLKMKGNFTLPEVAECFDEITYVELQKEEAQKLLEQYKEESKKALPPEKKQNTGSK KSNKNKSGKNQFNRRGGGHRGRGGFNMRGGNFRGGAPGNRGGYNRRGNMPQRGGGGGGSGGIGYPYPRAPVFPGRG SYSNRGNYNRGGMPNRRGNYNQNFRGRGNRRGYKNQSQGYNQWQQGQFWGQKPWSQHYHQGY

靶点信息

研究背景	This gene encodes a member of a family of proteins that bind nucleic acids and function in the formation of ribonucleoprotein complexes in the nucleus with heterogeneous nuclear RNA (hnRNA). The encoded protein has affinity for both RNA and DNA, and binds scaffold-attached region (SAR) DNA. Mutations in this gene have been associated with epileptic encephalopathy, early infantile, 54. A pseudogene of this gene has been identified on chromosome 14.
------	---

基因ID	3192
基因名	HNRNPU
Swiss	Q00839
别名	HNRNPU;EIEE54;HNRNPU-AS1;HNRPU;SAF-A;SAFA;U21.1;hnRNP U;hnRNP U;pp120

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn