

EEA1 Rabbit mAb

货号: B31079

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC
推荐浓度	WB: 1:500 - 1:2000 ICC: 1:50 - 1:200
理论分子量	162kDa
实测分子量	170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HT-1080,LO2,Mouse kidney,Mouse lung,Rat brain
细胞定位	Cytoplasm,Early endosome membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human EEA1 .
序列	LKAAVEQEKRNQILKDQVKKEEEEELKKEFIEKEAKLHSEIKEKEVGMKKHEENEAKLTMQITALNENLGTVKKEWQSSQRR VSELEKQTDDLERGEIAVLEATVQNNQDERRALLERCLKGEGEIEKLQTKVLELQRKLDNTTAAVQELGRENQSLQIKHTQAL NRKWAEDNEVQNCMACGKGFSVTVRRHHCRQCGNIFCAECSAKNALTPSSKKPVRVCDACFNDLQG

靶点信息

研究背景	Binds phospholipid vesicles containing phosphatidylinositol 3-phosphate and participates in endosomal tr afficking.
基因ID	8411

基因名	EEA1
Swiss	Q15075
别名	EEA1;MST105;MSTP105;ZFYVE2

产品验证



Western blot analysis of EEA1 expressed in Mouse kidney,Rat brain,U-87MG using EEA1 Rabbit mAb at 1 :1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn