

OLFM4 Rabbit mAb

货号: B31085

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	57kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse spleen,Mouse stomach,Rat stomach
细胞定位	Mitochondrion,Secreted,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human OLFM4.
序列	AIRREIVALKTKLKECEASKDQNTPVVHPPPTPGSCGHGGVNVISKPSVVQLNWRGFSYLYGAWGRDYSPQHHPNKGLYW VAPLNTDGRLLLEYRLYNTLDDLLLYINARELRITYGQSGTAVYNNNMVYVNMVNTGNIARVNLTTNTIAVTQTLPNAAAYNN RFSYANVAWQDIDFAVDENGLWVIYSTEASTGNMVISKLNDTTLQVLNTWYTKQYKPSASNAFMVCGVLYATRTMNTRT EEIFYYYDTNTGKEGKLDIVMHKMQEKVQSINYNPFDQKLYVYNDGYLLNYDLSVLQKPQ

靶点信息

研究背景	This gene was originally cloned from human myeloblasts and found to be selectively expressed in inflam med colonic epithelium. This gene encodes a member of the olfactomedin family. The encoded protein is an antiapoptotic factor that promotes tumor growth and is an extracellular matrix glycoprotein that facilit ates cell adhesion.
------	--

基因ID	10562
基因名	OLFM4
Swiss	Q6UX06
别名	OLFM4;GC1;GW112;OLM4;OlfD;UNQ362;bA209J19.1;hGC-1;hOLFd

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn