

GLUT1 Rabbit mAb

货号: B31086

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Cell membrane,Melanosome,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GLUT1.
序列	RPAAIAVAGFSNWTSTNFIVGMCQYVEQLCGPYVFIIFTVLLVLFIFTYFKVPETKGRTFDEIASGFRQGGASQSDKTPEEL FHPLGADSQV

靶点信息

研究背景	This gene encodes a major glucose transporter in the mammalian blood-brain barrier. The encoded protein is found primarily in the cell membrane and on the cell surface, where it can also function as a receptor for human T-cell leukemia virus (HTLV) I and II. Mutations in this gene have been found in a family with paroxysmal exertion-induced dyskinesia.
------	--

基因ID	6513
基因名	SLC2A1
Swiss	P11166
别名	SLC2A1;CSE;DYT17;DYT18;DYT9;EIG12;GLUT;GLUT-1;GLUT1;GLUT1DS;HTLVR;PED;SDCHCN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn