

HEC1 Rabbit mAb

货号: B31087

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	73kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,Lovo,MCF7,22Rv1,THP-1,293T,Mouse testis
细胞定位	Chromosome,Nucleus,centromere,kinetochore
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human HEC1.
序列	ETIKQENTRLQNIIDNQKYSVADIERINHERNELQQTINKLTKDLEAEQQKLWNEELKYARGKEAIETQLAEYHKLARKLKLIP KGAENSKGYDFEIKFNPEAGANCLVKYRAQYVPLKELLNETEEEINKALNKKMGLEDTLLEQLNAMITESKRSVRTLKEEVQ KLDDLYQQKIKEAEEEDKCASELESLEKHKHLLSTVQNQGLSEAMNELDAVQREYQLVVQTTTEERRKVGNNLQRLLEMV ATHVGSVEKHLEEQIAKVDREYEECMSEDLSENIKEIRDKYEKKATLIKSSEE

靶点信息

研究背景	This gene encodes a component of the NDC80 kinetochore complex. The encoded protein consists of an N-terminal microtubule binding domain and a C-terminal coiled-coiled domain that interacts with other components of the complex. This protein functions to organize and stabilize microtubule-kinetochore interactions and is required for proper chromosome segregation.
------	--

基因ID	10403
基因名	NDC80
Swiss	O14777
别名	NDC80;HEC;HEC1;HsHec1;KNTC2;TID3;hsNDC80

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn