

NUDEL Rabbit mAb

货号: B31088

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	37kDa/38kDa/40kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse brain
细胞定位	Chromosome,Cytoplasm,centromere,centrosome,cytoskeleton,kinetochore,microtubule organizing center ,spindle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NUDEL.
序列	MDGEDIPDFSSLKEETAYWKELSLKYQSFQEARDELVEFQEGSRELEAELEAQLVQAEQRNRDLQADNQRLKYEVEALK EKLEHQYAQSYKQVSVLEDDLSQTRAIKEQLHKYVRELEQANDDLERAKRATIVSLEDFEQRLNQAIERNAFLESELDEKES LLVSVQRLKDEARDLRQELAVRERQEQEVTRKSAPSSPTLDCEKMDSAVQASLSLPATPVGKGKTENTFPSPKAIPNGFGTSP LTPSARISALNIVGDLLRKVGALESKLAACRNFAKDQASRKSYSISGNVNCGLNNGNKGTKFSRSGHTSFFDKGAVNGFDPA PPPPGLGSSRPSSAPGMLPLSV

靶点信息

研究背景	This gene encodes a coiled-coil protein that plays a role in multiple processes including cytoskeletal orga nization, cell signaling and neuron migration, outgrowth and maintenance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is l ocated on the long arm of chromosome X.
------	---

基因ID	81565
基因名	NDEL1
Swiss	Q9GZM8
别名	NDEL1;EOPA;MITAP1;NDE1L1;NDE2;NUDEL

产品验证



Western blot analysis of NUDEL expressed in U-87MG,Mouse brain using NUDEL Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn