

PERP Rabbit mAb

货号: B31098

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	21kDa
实测分子量	21kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,K562,Jurkat,Mouse thymus
细胞定位	Cell junction,Cell membrane,Multi-pass membrane protein,desmosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PERP.
序列	CGPQMLVFLRVIGLLALAAVFQIISLVIYPVKYTQTFTLHANPAVTYIYNWAYGFGWAATIILIGCAFFFCCLPNYEDDLLGN AKPRYFYTSA

靶点信息

研究背景	Component of intercellular desmosome junctions. Plays a role in stratified epithelial integrity and cell-cell adhesion by promoting desmosome assembly. Plays a role as an effector in the TP53-dependent apoptotic pathway (By similarity).
基因ID	64065

基因名	PERP
Swiss	Q96FX8
别名	PERP;KCP1;KRTCAP1;PIGPC1;THW;dj496H19.1

产品验证



Western blot analysis of PERP expressed in Mouse thymus,HT-29,K562 using PERP Rabbit mAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn