

NDUFC2 Rabbit mAb

货号: B31123

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	10kDa/11kDa/13kDa/14kDa
实测分子量	12kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,OVCAR3,HeLa,A-549,mouse stomach,mouse brain,mouse heart,rat kidney
细胞定位	Matrix side,Mitochondrion inner membrane,Single-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NDUFC2.
序列	MIARRNPEPLRFLPDEARSLPPPKLTDPRLLYIGFLGYCSGLIDNLIIRRRPIATAGLHRQLLYITAFFFAGYYLVKREDYLYAVRD REMFGYMKLHPEDF

靶点信息

研究背景	Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I, th at is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH t o the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.
基因ID	4718

基因名	NDUFC2
Swiss	O95298
别名	NDUFC2;B14.5b;CI-B14.5b;HLC-1;NADHDH2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn