

RAP1A + RAP1B Rabbit mAb

货号: B31151

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP FC
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	21kDa
实测分子量	21kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,NIH/3T3,C6,U-251MG,HUVEC,Rat lung,Rat spleen
细胞定位	cell-cell junction,cytosol,extracellular exosome,glutamatergic synapse,lipid droplet,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human RAP1A + RAP1B.
序列	LYMKNGQG FALVYSITAQSTFNDLQDLREQILRVKDTDDVPMILVGKCDLEDERVVGKEQGQNLARQWNNCAFLESSA KS

靶点信息

研究背景	This gene encodes a member of the RAS-like small GTP-binding protein superfamily. Members of this family regulate multiple cellular processes including cell adhesion and growth and differentiation. This protein localizes to cellular membranes and has been shown to regulate integrin-mediated cell signaling. This protein also plays a role in regulating outside-in signaling in platelets. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 3, 5, 6 and 9. [provided by RefSeq, Oct 2011]
------	---

基因ID	5908,5906
基因名	RAP1B,RAP1A
Swiss	P61224/P62834
别名	K-REV; RAL1B

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn