

CHX10 Rabbit mAb

货号: **B31161**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	WB: 1:500 - 1:2000 FC: 1:20 - 1:50
理论分子量	50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spinal cord,Mouse brain,Rat spinal cord,Rat brain
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CHX10.
序列	AAAESGRKPEGERQALPKLDKMEQDERGPDAQAAISQEELRENSIAVLRKAQEHSTKVLGTVSGPDSLARSTEKPEEEEE MDEDRPAERLSPQLEDMA

靶点信息

研究背景	This gene encodes a homeobox protein originally described as a retina-specific transcription factor. Mutations in this gene are associated with microphthalmia, cataracts and iris abnormalities. [provided by RefSeq, Oct 2009]
基因ID	338917

基因名	VSX2
Swiss	P58304
别名	CHX10; HOX10; MCOP2; MCOPCB3; RET1

产品验证



Western blot analysis of CHX10 expressed in Mouse spinal cord,Rat spinal cord,Mouse brain using CHX10 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn