

PDE4B Rabbit mAb

货号: B31187

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	WB: 1:500 - 1:2000 FC: 1:20 - 1:50
理论分子量	57kDa/64kDa/82kDa/83kDa
实测分子量	95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,HL-60,Mouse brain,Mouse lung
细胞定位	centrosome,cytosol,dendritic spine,excitatory synapse,gamma-tubulin complex,nucleus,perinuclear region of cytoplasm,postsynaptic density,synaptic vesicle,Z disc
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PDE4B.
序列	MKKSRSVMTVMADDNVKDYFECSLSKSYSSSSNTLGIDLWRGRRCCSGNLQLPPLSQRQSERARTPEGDGISRPTTLPLT TLPSIAITTVSQECFDVENGPSPGRSPLDPQASSSAGLVLHATFPGHSQR

靶点信息

研究背景	This gene is a member of the type IV, cyclic AMP (cAMP)-specific, cyclic nucleotide phosphodiesterase (PDE) family. The encoded protein regulates the cellular concentrations of cyclic nucleotides and thereby plays a role in signal transduction. Altered activity of this protein has been associated with schizophrenia and bipolar affective disorder. Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding different isoforms.
------	---

基因ID	5142
基因名	PDE4B
Swiss	Q07343
别名	PDE4B;DPDE4;PDEIVB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn