

# VWF Rabbit mAb

货号: B31216

## 产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC                                                                                        |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200                                         |
| 理论分子量 | 309kDa                                                                                        |
| 实测分子量 | 309kDa                                                                                        |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | Mouse liver                                                                                   |
| 细胞定位  | Secreted,extracellular matrix,extracellular space                                             |
| 纯化    | Affinity purification                                                                         |

## 抗原信息


|      |                                                                                                                                                            |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human VWF.                                                                                                     |
| 序列   | APACQVTCQHPEPLACPVQCVEGCHAHCPPGKILDELLQTCVDPEDCPVCEVAGRRFASGKKVTLNPSDPEHCQICHCDV<br>VNLTCEACQEPGGLVVPPTDAPVSPPTLYVEDISEPPLHDFYCSRLDLVFLLDGSSRLSEAEFEVLKAFV |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a glycoprotein involved in hemostasis. The encoded preproprotein is proteolytically processed following assembly into large multimeric complexes. These complexes function in the adhesion of platelets to sites of vascular injury and the transport of various proteins in the blood. Mutations in this gene result in von Willebrand disease, an inherited bleeding disorder. An unprocessed pseudogene has been found on chromosome 22. [provided by RefSeq, Oct 2015] |
| 基因ID | 7450                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |

|       |            |
|-------|------------|
| 基因名   | VWF        |
| Swiss | P04275     |
| 别名    | F8VWF; VWD |

产品验证



Western blot analysis of VWF expressed in Mouse liver using VWF Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)