

Apg3 (Atg3) Rabbit mAb

货号: B31236

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	35kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,U-87MG,Mouse brain,Mouse kidney,Rat testis
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Apg3 (Atg3).
序列	MQNVINTVKGKALEVAEYLTPVLKESKFKETGVITPEEFVAAGDHLVHHCTWQWATGEELKVKAYLPTGKQFLVTKNVPCYKRCKQMEYSDELEAIIIEDDGGWVDYHNTGITGITEAVKEITLENKDNIRLQDCSALCEEEDEDEGEAADMEEYEE SGLLETDEATLDRKIVEACKAKTDAGGEDAILQTRTYDLYITYDKYYQTPRLWLFGYDEQRQPLTVEHMYEDISQDHVKK TVTIENHPLPPPPMCSVHPCRHAEVMKKIETVAEGGGELGVHMYLLIFLKFVQAVIPTIEYDYTRHFTM

靶点信息

研究背景	This gene encodes a ubiquitin-like-conjugating enzyme and is a component of ubiquitination-like systems involved in autophagy, the process of degradation, turnover and recycling of cytoplasmic constituents in eukaryotic cells. This protein is known to play a role in regulation of autophagy during cell death. A pseudogene of this gene is located on chromosome 20. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	--

基因ID	64422
基因名	ATG3
Swiss	Q9NT62
别名	ATG3;APG3;APG3-LIKE;APG3L;PC3-96

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn