

OS9 Rabbit mAb

货号: B31238

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	60kDa/63kDa/65kDa/67kDa/69kDa/73kDa/75kDa
实测分子量	83,97kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,LO2,BxPC-3,B cells,HeLa,Mouse liver
细胞定位	Endoplasmic reticulum lumen
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human OS9.
序列	NVQVKVIRSPADLIRFIEELKGGTKKGKPNIGQEQPVDDAAEVPQREPEKERGDPERQREMEEEDEDEDEDEDERQLL GEFEKELEGILLPSDRDRLRSEVKAGMERELENIIQETEKELDPDGLKKESERDRAMLALTSTLNKLIKRLLEEKQSPELVKKHK KKRVVPKKPPSPQPTGKIEIKIVRPWAEGTEEGARWLTDEDTRNLKEIFFNILVPGAEAAQKERQRQKELESNYRRVWGSP GGEGTGDLDEFDF

靶点信息

研究背景	This gene encodes a protein that is highly expressed in osteosarcomas. This protein binds to the hypoxia-inducible factor 1 (HIF-1), a key regulator of the hypoxic response and angiogenesis, and promotes the degradation of one of its subunits. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
------	--

基因ID	10956
基因名	OS9
Swiss	Q13438
别名	OS9;ERLEC2;OS-9

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn