

FOXF1 Rabbit mAb

货号: B31259

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	40kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse stomach
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human FOXF1.
序列	GLGFNHLPDTYGFQGSAGGLSCPPNSLALEGGGLGMMNGHLPGNVDGMALPSHSVPHLPNSNGGHSYMGCGGAAAGE YPHHDSSVPASPLLPTGAGGVMEPHAVYSGSAAAWPPSASAALNSGASYIKQQPLSPCNPAANPLSGSLSTHSLEQPYLH QNSHNAPAEHQGIPRYHSQSPSMCDRKEFVFSFNAMASSSMHSAGGGSYYHQVQTYQDIKPCVM

靶点信息

研究背景	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct fork head domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development.
基因ID	2294

基因名	FOXF1
Swiss	Q12946
别名	FOXF1;ACDMPV;FKHL5;FREAC1

产品验证



Western blot analysis of FOXF1 expressed in Mouse stomach using FOXF1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn