

IRS2 Rabbit mAb

货号: B31275

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	137kDa
实测分子量	185kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,K-562,Mouse kidney,Rat kidney
细胞定位	Cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human IRS2.
序列	LRAFCSHRSNTPESIAETPPARDGGGGGGEFYGYMTMDRPLSHCGRSYRRVSGDAAQDLDRGLRKRTYSLTTPARQRPVP QPSSASLDEYTLMRATFGSGAGRLCPSPASSPKVAYHPYPEDYGDIEIGSHRSSSSNLGAD

靶点信息

研究背景	This gene encodes the insulin receptor substrate 2, a cytoplasmic signaling molecule that mediates effect s of insulin, insulin-like growth factor 1, and other cytokines by acting as a molecular adaptor between div erse receptor tyrosine kinases and downstream effectors. The product of this gene is phosphorylated by t he insulin receptor tyrosine kinase upon receptor stimulation, as well as by an interleukin 4 receptor-asso ciated kinase in response to IL4 treatment.
------	---

基因ID	8660
基因名	IRS2
Swiss	Q9Y4H2
别名	IRS2;IRS-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn