

Phospho-GATA3 (S308) Rabbit mAb

货号: B31283

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	47kDa/48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Phospho-GATA3 (S308).
序列	MEVTADQPRWVSHHHPAVLNGQHPDTHHPGLSHSYMDAAQYPLPEEVDVLFNIDGQGNGHVPPYYGNSVRATVQRYPPT HHGSQVCRPPLLHGSLPWLDGGKALGSHHTASPWNLSPFSKTSIHGSGPGLSVYPPASSSSLSGGHASPFLTFPPTPP KDVSPDPSLSTPGSAGSARQDEKECLKYQVPLPDSMKLESSHSRGSM TALGGASSSTHHPITTPPYVPEYSSGLFPPSSL LGGSP TGF GCKSRPKARSSTG

靶点信息

研究背景	This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.
------	---

基因ID	2625
基因名	GATA3
Swiss	P23771
别名	GATA3;HDR;HDS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn