

SFRP4 Rabbit mAb

货号: B31315

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	39kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,HeLa,HO-8910
细胞定位	Secreted
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human SFRP4.
序列	VRGAPCEAVRIPMCRHMPWNITRMPNHLHHSTQENAILAIEQYEELVDVNCSAVLRFFLCAMYAPICTLEFLHDPIKPCKSVCC QRARDDCEPLMKMYNHSWPESLACDELVPYDRGVCSPEAIVTDLPEDVKWIDITPDMMVQERPLDVDCKRLSPDRCKCK KVKPTLATYLSKNYSYVIHAKIKAVQRSGCNEVTTVVDVKEIFKSSSPIRRTQVPLITNSSCQCPHILPHQDVLIMCYEWRSR MMLLENCLVEKWRDQLSKRSIQWEERLQEQRRTVQDKKKTAGRTSRSNPPKPKGKPPAPKPASPKNKIKTRSAQKRTNP KRV

靶点信息

研究背景	Secreted frizzled-related protein 4 (SFRP4) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. The expression of SFRP4 in ventricular myocardium correlates with apoptosis related gene expression.
基因ID	6424
基因名	SFRP4
Swiss	Q6FHJ7
别名	SFRP4;FRP-4;FRPHE;PYL;sFRP-4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn