

Exportin 5 Rabbit mAb

货号: B31375

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	136kDa
实测分子量	136kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,U-251MG,NIH/3T3,HeLa,PC-12,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Exportin 5.
序列	MRLSQKWQVINQRSLLCGEDEAADENPESQEMLEEQLVRLTREVMDLITVCCVSKKGADHSSAPPADGDDEEMMATE VTPSAMAELTDLGKCLMKHEDVCTALLITAFNSLAWKDTLSCQRTTSQLCWPLLKQVLSGTLADAVTWLFTSVLKGLQM HGQHDGCMASLVHLAFQIYEALRPRYLEIRAVMEQIPEIQKDSLQDFCKLLNPSLQKVADKRRKDQFKRLIAGCIGKPLGE QFRKEVHIKNLPSLFKKTCPMLETEVLNDGGGLATIFEP

靶点信息

研究背景	This gene encodes a member of the karyopherin family that is required for the transport of small RNAs and double-stranded RNA-binding proteins from the nucleus to the cytoplasm. The encoded protein translocates cargo through the nuclear pore complex in a RanGTP-dependent process.
------	--

基因ID	57510
基因名	XPO5
Swiss	Q9HAV4
别名	XPO5;exp5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn