

REST / NRSF Rabbit mAb

货号: B31376

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa/36kDa/119kDa/121kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Raji,SH-SY5Y,Mouse heart,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human REST / NRSF.
序列	ALPNDMYDLHDLSKAELAAPQLIMLANVALTGEVNGSCCDYLVGEERQMAELMPVGDNNFSDSEEGEGLEESADIKGEP HGLENMELRSLELSVVEPQPVFASGA

靶点信息

研究背景	This gene encodes a transcriptional repressor that represses neuronal genes in non-neuronal tissues. It is a member of the Kruppel-type zinc finger transcription factor family. It represses transcription by binding a DNA sequence element called the neuron-restrictive silencer element. The protein is also found in undifferentiated neuronal progenitor cells and it is thought that this repressor may act as a master negative regulator of neurogenesis. Alternatively spliced transcript variants have been described
基因ID	5978

基因名	REST
Swiss	Q13127
别名	REST;NRSF;WT6;XBR

产品验证



Western blot analysis of REST / NRSF expressed in HeLa,Raji,SH-SY5Y using REST / NRSF Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug p er lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 1 20s.

实验步骤

访问官网浏览详情: www.ablybio.cn