

GluR2+GluR3 Rabbit mAb

货号: B31390

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	98kDa
实测分子量	98kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	asymmetric synapse,dendrite,dendritic spine,excitatory synapse,external side of plasma membrane,plasma membrane,postsynapse,postsynaptic density,postsynaptic density membrane,postsynaptic endocytic zone
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GluR2+GluR3.
序列	KNKWWYDKGECGSGGGDSKEKTSALSLSNVAGVFYILVGGLGLAMLVALIEFCYSRAEAKRMKVAKNAQNINPSSSQNSQNFATYKEGYNVYGIESVKI

靶点信息

研究背景	
基因ID	2891,2892
基因名	GRIA2,GRIA3

Swiss	P42262/P42263
别名	

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn