

ATRIP Rabbit mAb

货号: B31402

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	72kDa/76kDa/83kDa/85kDa
实测分子量	86kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Jurkat,HeLa,Mouse brain,Mouse heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ATRIP.
序列	MAGTSAPGSKRRSEPPAPRPGPPPGTGHPPSKRARGFSAAAAPDPDDPFGAHGDFtADDLEELDTLASQALSQCpAAAR DVSSDHKvHRLLDGMSKNPSGKNRETVPiKDNFELEVLQAQYKELKEKMKVMEEVLIKNGEiKILRDSLHQTESVLEEQR RSHFLLEQEKtQALSdKEKEfSKKLQSLQSELQFKDAEMNELRTKLQTserANKLAAPSVSHVSPRKNPSVVIKPEACSPQ FGKTSFPTKESFSANMSLP

靶点信息

研究背景	This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangiecia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	84126
基因名	ATRIP
Swiss	Q8WXE1
别名	ATRIP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn