

TMOD2 Rabbit mAb

货号: B31403

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	35kDa/39kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse testis,Rat brain,Rat testis
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human TMOD2.
序列	EEALASASDTELYDLAAVLGVHNLNNPKFDEETANNKGGKGPVRNVVKGEKVKPVFEPPNPTNVEISLQQMKANDPSLQ

靶点信息

研究背景	This gene encodes a neuronal-specific member of the tropomodulin family of actin-regulatory proteins. The encoded protein caps the pointed end of actin filaments preventing both elongation and depolymerization. The capping activity of this protein is dependent on its association with tropomyosin. Alternatively spliced transcript variants encoding different isoforms have been described.
------	--

基因ID	29767
基因名	TMOD2
Swiss	Q9NZR1
别名	TMOD2;N-TMOD;NTMOD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn