

ADAR1 Rabbit mAb

货号: B31428

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	150kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T
细胞定位	Cytoplasm,Nucleus,Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ADAR1.
序列	STQAWNQHSGVVVRPDGHSQGAPNSDPSLEPEDRNSTSVSEDLLEPFIAVSAQAWNQHSGVVVRPDSHSQGSPNSDPGL EPEDSNSTSALEDPLEFLDMAEIK

靶点信息

研究背景	This gene encodes the enzyme responsible for RNA editing by site-specific deamination of adenosines. Th is enzyme destabilizes double-stranded RNA through conversion of adenosine to inosine. Mutations in this gene have been associated with dyschromatosis symmetrica hereditaria. Alternative splicing results in m ultiple transcript variants. [provided by RefSeq, Jul 2010]
------	--

基因ID	103
基因名	ADAR
Swiss	P55265
别名	ADAR1; AGS6; DRADA; DSH; DSRAD; G1P1; IFI-4; IFI4; K88DSRBP; P136

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn