

MTAP Rabbit mAb

货号: B31430

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	26kDa/30kDa/31kDa/32kDa/33kDa/36kDa/38kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29,A-549,NIH/3T3
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human MTAP.
序列	MASGTTTTAVKIGIIGGTGLDDPEILEGRTEKYVDTPFGKPSDALILGKIKNVDCVLLARHGRQHTIMPSKVNYQANIWALKE EGCTHVIVTTACGSLREEIQPGDIVIIDQFIDRTTMRPQSFYDGS HSCARGVCHIPMAEPFCPKTREVLIETAKKLGLRCHSK GTMVTIEGPRFSSRAESFMFRTWGADVINMTTVPEVVLAKEAGICYASIAMATDYDCWKEHEEAVSVDRVLKTLKENANK AKSLLLTTPIQIGSTEWSETLHNLKNMAQFSVLLPRH

靶点信息

研究背景	This gene encodes an enzyme that plays a major role in polyamine metabolism and is important for the salvage of both adenine and methionine. The encoded enzyme is deficient in many cancers because this gene and the tumor suppressor p16 gene are co-deleted. Multiple alternatively spliced transcript variants have been described for this gene, but their full-length natures remain unknown.
基因ID	4507
基因名	MTAP
Swiss	Q13126
别名	MTAP;BDMF;DMSFH;DMSMFH;HEL-249;LGMBF;MSAP;c86fus

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn